

Purge and Trap Troubleshooting: No Response, Carryover, and Contamination

Amy Nutter

VOC Technical Product Specialist



TELEDYNE TEKMAR
Everywhereyoulook™


Outline of Topics

- Overview of Purge and Trap analysis
- Troubleshooting
 - No Response
 - Carryover
 - Contamination




Teledyne Tekmar


- Created first commercial Purge and Trap (P&T) in 1975
- P&T Concentrators
 - Atomx XYZ
 - Lumin
 - AQUATek LVA




Brief Overview of Purge and Trap



Brief overview of Purge and Trap - Purge Mode



Brief Overview of Purge and Trap - Desorb Mode



Troubleshooting Tips

- Take a moment to review manuals
- Examine the facts and use valid reasoning
- Identify the “root cause” of the problem
- Avoid quick fixes and shotgun approach

BE PATIENT!



Troubleshooting – Finding the Root Cause

- MS, GC, Atomx XYZ (Concentrator and Autosampler)
 - Split the system into these 4 separate entities



Eliminate the GCMS

- Not necessary to disconnect P&T
- Confirm normal operation of the MS
 - Tune parameters/vacuum
 - Clean source
- Confirm normal operation of GC
 - Do direct inject of same standard used in P&T
- Poor chromatography/resolution/reduced response
 - Injection liner/septum
 - Clip column
- Once GCMS is eliminated as source of issue, proceed to P&T




Eliminate the P&T

- Confirm normal operation of the P&T
 - Run leak check
 - Make sure the correct analytical trap is in place
 - Manually load sample into P&T
 - Eliminates the autosampler




No Response

- Are the analytes getting to the GC/MS
 - Check the trap heater
 - Reaches the set temperature
 - Check 6-port valve
 - Hear the 6-port actuate
 - 6-port rotor installed correctly



No Response

- Are the analytes getting to the trap?
 - Is there a sample in the sparger and is it bubbling during purge?
 - Correct flow set in method
 - Mass flow controller working
 - System leak tight
 - Gas tank still have pressure



Carryover/Contamination

- Carryover for full target list/late eluters and contamination
- Run GC only
- Desorb only
- Check temperature
 - Desorb/Bake
- Check hot water/MeOH rinse



Carryover of Full Target List

- GC/MS
 - Faulty EPC
- P&T
 - Solenoid valve leaking cross port
 - Faulty trap heater
 - Sample not draining entirely
 - Dirty glassware
- Autosampler
 - Check soil vs. water
 - System not rinsing correctly
 - Internal standard valve leaking(allows standard to enter vessel)
 - Not transferring entire aliquot to P&T



Carryover of Late Eluters

- GC
 - Inlet too cold
- P&T
 - Not enough bake
 - Bad trap
 - Cold spots
 - Transfer line connection to GC inlet
 - Faulty condensate or trap heater
- Autosampler
 - Check soils vs. waters
 - Hot water heater/rinse failure
 - Not enough bake rinses
 - Sample needle not cleaning up
 - Not transferring bake rinse to P&T



How the Atomx XYZ Removes Contamination

- Able to rinse sample pathway with 90°C water
- Atomx XYZ uses methanol to rinse entire pathway
 - Includes syringe, associated transfer lines and valves, sparger and 3-stage needle.
 - User can vary volume and number of rinse times user settable within the software method



Thank You!

For more information:

Website: www.teledynetekmar.com
Phone: 513-229-7017 or 513-896-0702
Email: Amy.Nutter@Teledyne.com

Check out our website for all new applications!
You can also find us on Twitter, Facebook, and LinkedIn



TELEDYNE TEKMAR
Everywhereyoulook™