

Purge and Trap Workflows: Tips to Overcome Challenges That Can Arise From Sample Collection to Sample Analysis

Amy Nutter
Teledyne Tekmar
VOC Applications Chemist




TELEDYNE TEKMAR
Everywhereyoulook™

Outline of Topics


- Overview of Purge & Trap Technique and Tekmar Products
- Common Method Types and their Challenges
- How to Select Best Equipment for Method
- Consumables Considerations
- Tying it all Together



Brief overview of Purge and Trap




Brief overview of Purge and Trap





Introduction of Product Line

- Purge and Trap Concentrators
 - Atomx XYZ
 - Lumin
 - AQUATek LVA



Basic Method Types



Regulated Methods

- Globally, P&T usage is largely driven by regulated methods, all of which present unique challenges to end users
 - Sampling techniques and container considerations
 - Often, instrument conditions tightly controlled
 - Purge rates
 - Tuning requirements



What are the Challenges?

■ Vial cleanliness

- Contamination from vials can be frustrating to identify
 - All septa contribute siloxane peaks, worse in heated, in-vial purging
 - Target compounds are of a greater concern, can lead to failure of calibration and/or false positives
- Solution
 - Clean vials with soap and water, then bake before using
 - **Purchase certified pre-cleaned vials**
 - Keep sample heating to the minimum needed for low-level sensitivity



What are the Challenges?

- Cold preservation of samples
 - Once collected, must be kept on ice – makes shipments heavy
 - Once received into lab, must be kept in refrigerator until time of analysis
 - Some methods (EPA 524.3) require sample to be chilled on the autosampler as well.



What are the Challenges?

- Regulated parameters may not be the most efficient
 - USEPA (and China) method parameters are notorious for long run times and introducing significant amounts of water to GC
 - USEPA also regulates MSD parameters, often disallowing improvements in sensitivity by requiring BFB tune, not allowing SIM



What are the Challenges?

EPA 524.2

Parameter	Time (min)
Purge	11
desorb	4
Bake	4
Total	19!

EPA 8260


Parameter	Time (min)
Purge	11
desorb	2
Bake	4
Total	17!

EPA 524.3


Parameter	Time (min)
Purge	5 to 11
desorb	0.5 to 2
Bake	2
Total	7.5 to 15!

What are the Challenges?

- Methods introduce large amounts of water to GC



Water Reduction per Desorb Time



Desorb Time (min)	% of 4 min
4	100
3	93
2	58
1	28
0.5	13

Method-specific Challenges

- Each method, or group of methods has their own collection of challenges.
 - Drinking water – low level sensitivity
 - Wastewater – carryover and dynamic range
 - Specialty methods
 - Mold/Odor analysis – purge efficiency
 - 1,4-dioxane – low level sensitivity



Drinking Water Method Challenges

- Low level sensitivity
 - As source, column, trap and sparger become dirty, low-end sensitivity becomes an issue
 - Moisture control is particularly important to controlling this
 - Cleanliness of purge gas, rinse water and glassware associated with standards prep
 - Sparger and trap maintenance



Wastewater Method Challenges

- Carryover can cause false positives or require additional rinses or blanks – increase run time
- Source, column, trap and sparger maintenance
- Active sites in gas pathway
 - Activity in tubing, trap and MCS can all contribute to poor calibrations and loss of compounds
 - Clean and/or replace pathway



Specialty Methods Challenges

■ Mold/Odor

■ Method parameters often require optimization

- Application note is a good starting point, not always workable solution
- Can benefit from use of a sparge vessel heater

■ Compounds are large molecules, tend to foul traps and cause backpressure issues

- Recommend #1A trap, large mesh tenax offers decent sensitivity while lasting longer than other options



Specialty Methods Challenges

- 1,4-dioxane
 - Not very volatile, shows poor purging characteristics
 - Often requires large sample volumes (25ml sparger)
 - SIM scanning can be helpful in improving sensitivity
 - Application can benefit from use of sparge vessel heater
 - K or #9 trap are sufficient for this analysis

What are the Solutions?

- Improvements in Tekmar hardware have allowed some improvements in time
 - Faster trap cool down from bake to purge ready
 - Faster trap heating shortens time to reach desorb temp.
 - Improvements in moisture management reduce or eliminate need for dry purge



What are the Solutions?

- Moisture management is controlled through hardware and application parameters
 - Improved moisture control fitting in Lumin and Atomx XYZ
 - Shorter, narrow bore GC columns allow for higher split ratios, reducing moisture introduced
 - Shorter desorb times reduce water introduced

Packages can be Broken into Key Methods

■ Drinking Water

- Sparger purge
 - Need chilling
 - Don't need chilling
- In-vial purge


■ Wastewater/soils

- Need in-vial purge capabilities
- Do not need in-vial purge
 - Need dilutions or MeOH extraction
 - With or without antifoam
 - No need for dilutions or MeOH extraction
 - With or without antifoam




Drinking Water Packages

- Lumin and AQUATek LVA
 - Offers chiller standard, 84 vial positions, fixed sample loop (available in multiple sizes), same concentrator tech found in Atomx XYZ
 - K (vocarb 3000) or #9 trap




Drinking Water Packages

- Atomx XYZ
 - Typically, these users are not interested in dilutions, MeOH rinse or extractions
 - K (vocarb 3000) or #9 trap




Wastewater Packages

- Atomx XYZ
 - Offers dilutions, soil purge, MeOH extraction and rinsing all in base model
 - Top model add foam eliminator and chiller
 - K (vocarb 3000) or #9 trap
 - Consider adding **fritless sparger** if foam is a concern



Wastewater Packages

- Sparger Purge
 - Lumin and AQUATek LVA
 - Good for customers who do not analyze soils, and do not need dilutions capability.
 - If foam is a concern, foam eliminator model of Lumin and/or **fritless sparger**
 - K (vocarb 3000) or #9 trap



Thank You!

For more information:

Website: www.teledynetekmar.com

Phone: 513-229-7017

Email: amy.nutter@Teledyne.com

Check out our website for all new applications!

You can also find us on Twitter, Facebook, and LinkedIn



TELEDYNE TEKMAR
Everywhereyoulook™