

Troubleshooting Any Purge and Trap

Amy Nutter
VOC Technical Product Specialist
March 2023



TELEDYNE TEKMAR
Everywhereyoulook™

Outline of Topics


- Overview of Teledyne Tekmar and Our Products
- Troubleshooting
- TekLink Troubleshooting Tools
- Questions?



TELEDYNE TEKMAR
Everywhereyoulook™

Teledyne Tekmar

- Created first commercial Purge and Trap (P&T) in 1975
- P&T Concentrators
 - Atomx XYZ
 - Lumin
 - AQUATek LVA



P&T Troubleshooting

- No Response
- Carryover/Contamination
- Elution issue- early, middle, late
- Water interference
- Internal Standard recovery variability
- Linearity/Reproducibility issues
- System leak check
- Benchmark test



Troubleshooting Tips

- Take a moment to review manuals
- Examine the facts and use valid reasoning
- Identify the “root cause” of the problem
- Avoid quick fixes and shotgun approach

BE PATIENT!



TELEDYNE TEKMAR
Everywhereyoulook™

Troubleshooting – Finding the Root Cause

- MS, GC, Concentrator, and Autosampler
 - Split the system into these 4 separate entities



Eliminate the GCMS

- Not necessary to disconnect P&T
- Confirm normal operation of the MS
 - Tune parameters/vacuum
 - Clean source
- Confirm normal operation of GC
 - Do direct inject of same standard used in P&T
- Poor chromatography/resolution/reduced response
 - Injection liner/septum
 - Clip column
- Once GCMS is eliminated as source of issue, proceed to P&T



Eliminate the P&T

- Confirm normal operation of the P&T
 - Run leak check
 - Make sure the correct analytical trap is in place
 - Manually load sample into P&T
 - Eliminates the autosampler



No Response

- Are the analytes getting to the GC/MS
 - Check the trap heater
 - Reaches the set temperature
 - Check 6-port valve
 - Hear the 6-port actuate
 - 6-port rotor installed correctly



TELEDYNE TEKMAR
Everywhereyoulook™

No Response

- Are the analytes getting to the trap?
 - Is there a sample in the sparger and is it bubbling during purge?
 - Correct flow set in method
 - Mass flow controller working
 - System leak tight
 - Gas tank still have pressure



Carryover/Contamination

- Carryover for full target list/late eluters and contamination
- Run GC only
- Desorb only
- Check temperature
 - Desorb/Bake
- Check hot water/MeOH rinse



Carryover of Full Target List

- GC/MS

- Faulty EPC

- P&T

- Solenoid valve leaking cross port
 - Faulty trap heater
 - Sample not draining entirely
 - Dirty glassware

- Autosampler

- Check soil vs. water
 - System not rinsing correctly
 - Internal standard valve leaking(allows standard to enter vessel)
 - Not transferring entire aliquot to P&T



Carryover of Late Eluters

- GC
 - Inlet too cold
- P&T
 - Not enough bake
 - Bad trap
 - Cold spots
 - Transfer line connection to GC inlet
 - Faulty condensate or trap heater
- Autosampler
 - Check soils vs. waters
 - Hot water heater/rinse failure
 - Not enough bake rinses
 - Sample needle not cleaning up
 - Not transferring bake rinse to P&T



Elution Issues

- Tune
 - Auto-tune/BFB tune
- Direct inject
 - If missing peaks could be a bad EPC
 - Column
- Top of Trap
 - If missing peaks could be sample tubing/transfer line



Early Elution Groups

- GC

- Split flow leak
- Bad column

- P&T

- Small leak
- Water management
- Too long of dry purge
- Trap too hot (faulty electronics or heater)



Middle Elution Groups (Very Rare)

- GC
 - Column phase is bad
- P&T
 - Desorb preheat problems
 - Mass flow controller leaking
 - Bad Trap



Late Elution Groups

- GCMS
 - Source maintenance
 - Carrier flow too low
- P&T
 - Poor purge efficiency
 - Restrictive sample tubing
 - Dirty glassware
 - Sample not heated properly
 - Plugged vent valve
- Bad Trap



Water Interference

- Purge and Trap
 - Trap type
 - Sample purge temperature
 - Dry purge settings not optimized
 - Desorb time
 - Bake time and temperature



Water Interference

- GC
 - Split ratio
 - Might be time to replace column
- MS
 - Vacuum issues



Internal Standard Recoveries

- GCMS response increases as Target compound concentration increases
 - Dirty inlet liner or source
 - Validate by preparing 3 increasing target concentrations with Internal Standard concentration fixed in 1mL vials and then direct inject 1uL of the concentrations into the GC
 - If still increasing isolate between source and liner by performing cleaning and or replacement.



Internal Standard Recoveries

- P&T
 - Isolated to an active site
 - Sample tubing/analytical trap
- Autosampler
 - Hand spike vials first
 - Automated delivery system
 - Tubing restrictions
 - Valve malfunction
 - Pressure problems



Linearity and Reproducibility

■ GCMS

- Source maintenance
- Vacuum problems
- EPC failures
- Bad column
- Method not optimized
 - Oven temperature program
- Multiplier going bad



Linearity and Reproducibility

- Purge and Trap
 - Method not optimized
 - Trap is failing
 - Drain valve leaking
 - Excess water in system
 - Temperatures not reaching set-points
 - Too much carry over
 - Sample preparation errors
 - Dirty standard syringe or glassware
 - Too much headspace
 - Analyst variability



Linearity and Reproducibility

■ Autosampler

- Not pulling same amount of sample
- Not transferring sample correctly
- Improper rinsing between samples
- Compare soil vs. water samples
- Determine if IS is varying and causing everything else to look unstable




Tips to Identify Leaks Before They Become a Problem: Warning Signs

- Check History Log in Tools and review purge and bake pressures
- Front end compounds are missing/low response
- Linearity and reproducibility
- Slow draining of the sparger
- Low/No sparge bubbles
- GC beeping/alarming during desorb




Tips to Identify Leaks Before They Become a Problem: Warning Signs




TekLink™ Leak Check Diagnostic Tool

- Perform a Leak check:
 - after changing the IS, foam eliminator, any glassware, trap, and/or any fittings
 - after turning the system back on after being off for a while
 - before performing a benchmark test



Benchmark

- Make sure system is leak tight first
- Checks all system controls including solenoid valves, heaters, flows, and 6-port valve
- Helps isolate a variety of potential issues



Thank You!

For more information:

Website: www.teledynetekmar.com
Phone: 513-229-7017 or 513-896-0702
Email: Amy.Nutter@Teledyne.com

Check out our website for all new applications!
You can also find us on Twitter, Facebook, and LinkedIn



TELEDYNE TEKMAR
Everywhere you look™